


medax

CAESAR

AUTOMATIC DISPOSABLE BIOPSY SYSTEM



Automatic disposable spring loaded biopsy system.

- One hand use.
- 10/22 mm penetration depth.
- Safety switch automatically activated.
- Double firing trigger in proximal and distal positions.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CC introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

AUTOMATIC DISPOSABLE BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)							
14	● Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
16	○ White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300
18	● Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	● Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300

Different sizes on demand

Ordering information:

CAESAR NORMAL TIP (penetration depth 22 mm)

C	S	GAUGE		LENGTH		0	0

CAESAR TROCAR TIP (penetration depth 22 mm)

C	S	GAUGE		LENGTH		0	1

CAESAR SHORT NOTCH NORMAL TIP (penetration depth 10 mm)

C	S	GAUGE		LENGTH		0	S

INTRODUCER COAXIAL NEEDLE FOR CAESAR

C	C	GAUGE		LENGTH		0	0

CAESAR WITH COAXIAL NEEDLE IN THE SAME PACKAGE

C	S	GAUGE		LENGTH		C	0

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

