

medax

MEDONE

AUTOMATIC DISPOSABLE BIOPSY SYSTEM



Automatic disposable spring loaded biopsy system.






- One hand usage.
- 22 mm penetration depth.
- Double firing trigger in proximal and distal positions.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CO introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

MEDONE

AUTOMATIC DISPOSABLE BIOPSY SYSTEM

| GAUGE | COLOUR | LENGTH (MM) | | | | | | | |
|-------|--|-------------|--------|--------|--------|--------|--------|--------|--------|
| 12 |  Light blue | 12 080 | 12 100 | 12 120 | 12 160 | 12 180 | 12 200 | 12 250 | 12 300 |
| 14 |  Green | 14 080 | 14 100 | 14 120 | 14 160 | 14 180 | 14 200 | 14 250 | 14 300 |
| 16 |  White | 16 080 | 16 100 | 16 120 | 16 160 | 16 180 | 16 200 | 16 250 | 16 300 |
| 18 |  Pink | 18 080 | 18 100 | 18 120 | 18 160 | 18 180 | 18 200 | 18 250 | 18 300 |
| 20 |  Yellow | 20 080 | 20 100 | 20 120 | 20 160 | 20 180 | 20 200 | 20 250 | 20 300 |

Different sizes on demand

Ordering information:

MEDONE NORMAL TIP

| | | | | | | | | |
|---|---|----------------------|----------------------|----------------------|----------------------|----------------------|---|---|
| M | O | GAUGE | | LENGTH | | | 0 | 0 |
| | | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | | |

MEDONE TROCAR TIP

| | | | | | | | | |
|---|---|----------------------|----------------------|----------------------|----------------------|----------------------|---|---|
| M | O | GAUGE | | LENGTH | | | 0 | 1 |
| | | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | | |

INTRODUCER COAXIAL NEEDLE FOR MEDONE

| | | | | | | | | |
|---|---|----------------------|----------------------|----------------------|----------------------|----------------------|---|---|
| C | O | GAUGE | | LENGTH | | | 0 | 0 |
| | | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | | |

MEDONE WITH COAXIAL NEEDLE IN THE SAME PACKAGE

| | | | | | | | | |
|---|---|----------------------|----------------------|----------------------|----------------------|----------------------|---|---|
| M | O | GAUGE | | LENGTH | | | C | O |
| | | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | | |

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

