


medax

MEDEXTRA

LYMPH NODE BIOPSY SYSTEM



medax
medical devices


Your partner in Biopsy and Special Needles

medax

MEDEXTRA

LYMPH NODE BIOPSY SYSTEM


- Patented rectangular shaped cannula.
- 4 mm wide flat biopsy sample.
- Ultrasharp stylet tip clearly visible with ultrasound guidance.
- Adjustable penetration depth 12/22 mm.
- Centimeter markings.
- Ergonomic handle.
- Easy sample retrieval and release.
- Sterilized by ETO, shelf life 5 years.




LYMPH NODE BIOPSY SYSTEM designed to generate large amount of sample with low bleeding risk and be well tolerated by patients, for quicker diagnosis at critical time and implementation of earlier treatment.

FEATURES


- 4 mm wide flat biopsy sample (4x20x1 mm).
- 3 times surface area and volume of sample compared to standard core biopsy needles.
- Plenty of cleanly cut sample for histopathologist and for further specialist analyses.
- Broad architectural view of tissue samples under the microscope.
- Special tip designed to help easy and a-traumatic penetration.




Regular core biopsy sample



Medextra wide flat sample
(4x20x1mm)



Highly visible under
ultrasound guidance



Sharp cutting edges stylet

BENEFITS

- Improved shape leads to bigger tissue quantity for special analyses and diagnosis.
- A-traumatic biopsy procedure helps patients to recover earlier.
- Suitable for use under local anaesthetic avoiding general anaesthesia risks.
- Lower cost compared to whole node excision.

Ordering information:

MEDEXTRA

X T 0 0 1 0 0 - 0 0

Single box: 5 PCS.

WHEN IT IS URGENT TO OBTAIN GOOD TISSUE IN SUSPECTED LYMPHOMA AND YOU WANT TO AVOID THE RISKS AND DELAYS IN ARRANGING A GENERAL ANAESTHETIC, THEN GO FOR MEDEXTRA BIOPSY SYSTEM.

DESIGNED BY DOCTORS WHO UNDERSTAND DOCTORS



Your partner in Biopsy and Special Needles

Medax Srl Unipersonale

Via S. Pertini, 4 • 41039 • San Possidonio (MO) Italy

Tel. +39 0535 1812757 • Fax +39 0535 1812744

www.medax.it • customercare@medax.it



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.