


medax

UNICORE

MEDCORE COMPATIBLE BIOPSY NEEDLE
ALSO USABLE WITH BARD® MAGNUM®
AUTOMATIC REUSABLE BIOPSY SYSTEM



- MEDCORE compatible needle.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CU introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

UNICORE

MEDCORE COMPATIBLE BIOPSY NEEDLE
 ALSO USABLE WITH BARD® MAGNUM® AUTOMATIC REUSABLE BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)							
12	● Light blue	12 080	12 100	12 120	12 160	12 180	12 200	12 250	12 300
14	● Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
16	○ White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300
18	● Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	● Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300

Different sizes on demand

Ordering information:

UNICORE NORMAL TIP

U	C	GAUGE		LENGTH			0	0

UNICORE TROCAR TIP

U	C	GAUGE		LENGTH			0	1

INTRODUCER COAXIAL NEEDLE FOR UNICORE

C	U	GAUGE		LENGTH			0	0

UNICORE WITH COAXIAL NEEDLE IN THE SAME PACKAGE

U	C	GAUGE		LENGTH			C	0

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

