


**medax**

# UNIVERSAL

M-GUN COMPATIBLE BIOPSY NEEDLE  
ALSO USABLE WITH MOST COMMON MARKET  
REUSABLE BIOPSY SYSTEMS



M-Gun compatible needle also usable with the most common reusable guns:





- Angiotech/Pro-mag™ 2.2
- Bard®/Biopty®
- Gallini/ABS
- US Biopsy/MBD-23
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CU introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

# UNIVERSAL

M-GUN COMPATIBLE BIOPSY NEEDLE ALSO USABLE WITH AND MOST COMMON MARKET REUSABLE BIOPSY SYSTEMS

GAUGE	COLOUR	LENGTH (MM)								
14	 Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300	
16	 White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300	
18	 Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300	
20	 Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300	

Different sizes on demand

Ordering information:

UNIVERSAL NORMAL TIP

U	N	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

UNIVERSAL TROCAR TIP

U	N	GAUGE		LENGTH			0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

INTRODUCER COAXIAL NEEDLE FOR UNIVERSAL

C	U	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

UNIVERSAL WITH COAXIAL NEEDLE IN THE SAME PACKAGE

U	N	GAUGE		LENGTH			C	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

