

medax

AMNIOX

SPECIAL SHARPENED NEEDLES FOR TRANSABDOMINAL
RETRIEVAL OF AMNIOTIC FLUID
DURING PRENATAL DIAGNOSIS



Special sharpened needles for transabdominal retrieval of amniotic fluid during prenatal diagnosis.

- Transparent Luer Lock connector.
- Sharpened cannula for a safe and painless insertion and penetration.
- Internal or external echogenic treatment for a correct positioning under ultrasound guidance.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

MEDAX

AMNIOX

SPECIAL SHARPENED NEEDLES FOR TRANSABDOMINAL RETRIEVAL OF AMNIOTIC FLUID DURING PRENATAL DIAGNOSIS

GAUGE	COLOUR	LENGTH (MM)		
20	 Yellow	20 120	20 160	20 180
21	 Light green	21 120	21 160	21 180
22	 Black	22 120	22 160	22 180

Ordering information:

AMNIOX



Single box: 20 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.