


**medax**

**COX**

## STANDARD COAXIAL INTRODUCER NEEDLE



Standard coaxial introducer needle used for single access in multiple biopsies.

- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Standard Trocar tip.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

**MEDAX****COX****STANDARD COAXIAL INTRODUCER NEEDLE**

BIOPSI NEEDLE GAUGE	COX COLOUR	COX LENGTH (MM)							
		80	100	120	160	180	200	250	300
12	● Light green	11 030	11 050	11 070	11 110	11 130	11 150	11 200	11 250
14	● Light blue	13 030	13 050	13 070	13 110	13 130	13 150	13 200	13 250
16	● Blue	15 030	15 050	15 070	15 110	15 130	15 150	15 200	15 250
18	● Red	17 030	17 050	17 070	17 110	17 130	17 150	17 200	17 250
20	● Light brown	19 030	19 050	19 070	19 110	19 130	19 150	19 200	19 250

To find the desired code for the coaxial needle, cross the measures of the biopsy needle displayed in the black boxes (if you need a coaxial for a 16G100mm needle, cross the reference to find the correct coaxial code, in this case 15050).  
Different sizes on demand

**Ordering information:**

COX



Single box: 20 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.