

**medax**

**FNA**

## FINE ASPIRATION NEEDLES FOR CYTOLOGICAL AND MICROHISTOLOGICAL BIOPSY

**CHIBA**




**QUINCKE**




**TURNER**



**WESTCOTT**



**FRANSEEN**



Fine aspiration needles for cytological and microhistological biopsy with different tip versions: CHIBA, QUINCKE, WESTCOTT, FRANSEEN and TURNER.

- Transparent Luer Lock connector.
- Sharpened cannula for a safe and painless insertion and penetration.
- Internal or external echogenic treatment for a correct positioning under ultrasound guidance.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.

### **CHIBA and QUINCKE needle:**

Disposable needle used for cytological aspiration and local anaesthesia administration.

### **WESTCOTT, FRANSEEN and TURNER needle:**

Disposable needle used for cytological/microhistological biopsy.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

## FINE ASPIRATION NEEDLES FOR CYTOLOGICAL AND MICROHISTOLOGICAL BIOPSY

GAUGE	COLOUR	LENGTH (MM)							
18	Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300
21	Light green	21 080	21 100	21 120	21 160	21 180	21 200	21 250	21 300
22	Black	22 080	22 100	22 120	22 160	22 180	22 200	22 250	22 300

GAUGE	COLOUR	LENGTH (MM)						
23	Light blue	23 040	23 050	23 070	23 090	23 110	23 150	23 200
25	Orange	25 040	25 050	25 070	25 090	25 110	25 150	25 200

Different sizes on demand - WESTCOTT and FRANSEEN not available for Gauge 23 and Gauge 25.

### Ordering information:

CHIBA NEEDLE

C	H	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

CHIBA NEEDLE WITH ECHOGENIC TREATMENT

C	H	GAUGE		LENGTH			0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

QUINCKE NEEDLE

Q	K	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

QUINCKE NEEDLE WITH ECHOGENIC TREATMENT

Q	K	GAUGE		LENGTH			0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

WESTCOTT NEEDLE

W	E	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

WESTCOTT NEEDLE WITH ECHOGENIC TREATMENT

W	E	GAUGE		LENGTH			0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

FRANSEEN NEEDLE

F	R	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

FRANSEEN NEEDLE WITH ECHOGENIC TREATMENT

F	R	GAUGE		LENGTH			0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

TURNER NEEDLE

T	U	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

TURNER NEEDLE WITH ECHOGENIC TREATMENT

T	U	GAUGE		LENGTH			0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

Single box: 20 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

