


**medax**

**FACILE**

TRU-CUT MANUAL BIOPSY SYSTEM



Tru-cut manual biopsy system

- 22 mm penetration depth.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CF introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

**MEDAX**

# FACILE

## TRU-CUT MANUAL BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)							
		14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
14	● Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
16	○ White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300
18	● Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	● Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300
21	● Light green	21 080	21 100	21 120	21 160	21 180	21 200	21 250	21 300

Different sizes on demand

### Ordering information:

FACILE NORMAL TIP

F	C	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

FACILE TROCAR TIP

F	C	GAUGE		LENGTH			0	1
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

INTRODUCER COAXIAL NEEDLE FOR FACILE

C	F	GAUGE		LENGTH			0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

FACILE WITH COAXIAL NEEDLE IN THE SAME PACKAGE

F	C	GAUGE		LENGTH			C	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.