


medax

HANDLEX

BONE MARROW ASPIRATION SYSTEM
WITH ERGONOMIC HANDLE "HN"
ALSO AVAILABLE IN THE EXPLANT/TRANSPLANT
VERSION "HD"



EXPLANT/TRANSPLANT VERSION

Bone Marrow aspiration system with ergonomic handle "HN", also available in the explant/transplant version "HD".

- Ergonomic handle for a comfortable procedure and a firm grip.
- Triple sharpened tip cannula to allow a non traumatic, fast and easy penetration into the bone marrow cavity.
- Adjustable depth setter with sliding mechanism for fast and precise settings.
- Luer Lock connection for syringe.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

HANDLEX

BONE MARROW ASPIRATION SYSTEM WITH ERGONOMIC HANDLE "HN"
 ALSO AVAILABLE IN THE EXPLANT/TRANSPLANT VERSION "HD"

Ordering information:

| GAUGE | COLOUR | LENGTH (MM) | | | | | |
|-------|---------|-------------|--------|--------|--------|--------|--------|
| 14 | ● Green | 14 030 | 14 040 | 14 050 | 14 070 | 14 090 | 14 110 |
| 15 | ● Blue | 15 030 | 15 040 | 15 050 | 15 070 | 15 090 | 15 110 |
| 16 | ○ White | 16 030 | 16 040 | 16 050 | 16 070 | 16 090 | 16 110 |
| 18 | ● Pink | 18 030 | 18 040 | 18 050 | 18 070 | 18 090 | 18 110 |

Different sizes on demand

Ordering information:

HANDLEX ASPIRATION TIP

| | | | | | | | | |
|---|---|-------|--|--------|--|--|---|---|
| H | N | GAUGE | | LENGTH | | | 0 | 0 |
|---|---|-------|--|--------|--|--|---|---|

HANDLEX EXPLANT/TRANSPLANT TIP

| | | | | | | | | |
|---|---|-------|--|--------|--|--|---|---|
| H | D | GAUGE | | LENGTH | | | 0 | 0 |
|---|---|-------|--|--------|--|--|---|---|

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

