

medax

HEPAX

**HISTOLOGICAL BIOPSY
MODIFIED MENGHINI INSTRUMENT**



Ultrasound guided instrument for histological aspiration biopsy using the modified Menghini technique, recommended for a safe and painless liver biopsy.

- Syringe of 10 cc with self-locking system to guarantee a constant vacuum during the procedure.
- Ultra thin sharpened cannula with trocar stylet for a quick and painless percutaneous insertion.
- Internal (on the stylet) or external echogenic treatment for a correct positioning under ultrasound guidance.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

MEDaX

HEPAX

HISTOLOGICAL BIOPSY MODIFIED MENGHINI INSTRUMENT

GAUGE	COLOUR	LENGTH (MM)				
15	● Blue	15 080	15 100	15 120	15 160	15 200
16	○ White	16 080	16 100	16 120	16 160	16 200
17	● RED	17 080	17 100	17 120	17 160	17 200
18	● Pink	18 080	18 100	18 120	18 160	18 200
20	● Yellow	20 080	20 100	20 120	20 160	20 200

Different sizes on demand

Ordering information:

HEPAX



Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.