


medax

MRI KIT

SEMI-AUTOMATIC SPRING LOADED BIOPSY SYSTEM
WITH MRI COMPATIBLE INTRODUCER COAXIAL NEEDLE



MRI compatible coaxial
introducer needle

BIO-FEATHER

Semi-Automatic spring loaded
biopsy system

MRI KIT composed by a semi-automatic spring loaded biopsy system
and a MRI compatible introducer coaxial needle.

- Kit usable in MRI biopsy procedures.
- Centimeter markings and optional sliding stopper.
- Sterilized by ETO, shelf life 5 years.

MRI compatible introducer coaxial needle.

- Usable in MRI biopsy procedures.
- Centimeter markings and optional sliding stopper.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

MRI KIT

SEMI-AUTOMATIC SPRING LOADED BIOPSY SYSTEM
WITH MRI COMPATIBLE INTRODUCER COAXIAL NEEDLE

GAUGE	COLOUR	LENGTH (MM)		
14	● Green	14 100	14 160	14 200
16	○ White	16 100	16 160	16 200
18	● Pink	18 100	18 160	18 200

Different sizes on demand

Ordering information:

MRI KIT



Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.