


medax

MEDEM

MALLARME BONE MARROW ASPIRATION SYSTEM
WITH ERGONOMIC HANDLE "MD"
ALSO AVAILABLE IN THE EXPLANT/TRANSPLANT VERSION "MX"



Mallarmé Bone Marrow aspiration system with ergonomic handle "MD".
Also available in the explant/transplant version "MX"

- Ergonomic handle for a comfortable procedure and a firm grip.
- Triple sharpened tip cannula to allow a non traumatic, fast and easy penetration into the bone marrow cavity.
- Luer Lock connection for syringe.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

MEDEM

MALLARME BONE MARROW ASPIRATION SYSTEM WITH ERGONOMIC HANDLE "MD"
 ALSO AVAILABLE IN THE EXPLANT/TRANSPLANT VERSION "MX"

Ordering information:

GAUGE	COLOUR	LENGTH (MM)							
14	● Green	14 020	14 030	14 035	14 040	14 045	14 050	14 055	14 060
15	● Blue	15 020	15 030	15 035	15 040	15 045	15 050	15 055	15 060
16	○ White	16 020	16 030	16 035	16 040	16 045	16 050	16 055	16 060
18	● Pink	18 020	18 030	18 035	18 040	18 045	18 050	18 055	18 060

Different sizes on demand

Ordering information:

MEDEM

M	D	GAUGE		LENGTH			0	0
---	---	-------	--	--------	--	--	---	---

MEDEM EXPLANT/TRANSPLANT TIP

M	X	GAUGE		LENGTH			0	0
---	---	-------	--	--------	--	--	---	---

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

