

medax

MEDLOCK

BONE MARROW BIOPSY AND ASPIRATION SYSTEM
WITH SPECIAL EXTRACTION CANNULA



EXTRACTION CANNULA
"MF" VERSION

EXTRACTION CANNULA
"MM" VERSION



Bone marrow biopsy and aspiration system provided
with special safe and atraumatic extraction cannula.

- Coloured-code ergonomic handle.
- Additional extraction cannula to get the sample in a completely atraumatic and safe way.
- "MM" version optimal for multiple procedures.
- Luer Lock connection for syringe.
- Remover guide provided with indicators to facilitate the sample expulsion and to check its length during the procedure.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

MEDLOCK

**BONE MARROW BIOPSY AND ASPIRATION SYSTEM
WITH SPECIAL EXTRACTION CANNULA**

GAUGE	COLOUR	LENGTH (MM)		
8	● Dark Green	08 070	08 100	08 150
9	● Pink	09 070	09 100	09 150
11	○ White	11 070	11 100	11 150
13	● Light blue	13 070	13 100	13 150

Different sizes on demand

Ordering information:

MEDLOCK "MM" VERSION DIAMOND TIP

M	M	GAUGE	LENGTH	0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

MEDLOCK "MF" VERSION FISH MOUTH TIP

M	F	GAUGE	LENGTH	0	0
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.