

medax

TORASET

SET FOR THORACENTESIS



Toracentesis Set for aspiration of pleuric exudates.

- Practical system with single-acting valves.
- Waste bag with non-return valve and drainage tap.
- Available in various kit configurations.
- Wide range of needles:
 - Trio of 55 mm needles in the diameters 14 G, 16 G, 19 G;
 - Trio of 80 mm needles in the diameters 14 G, 16 G, 19 G;
 - Veress needle.
- Available in manual version (TORASET) with a tightly sealed 3-way stopcock on the connecting tubing, or in automatic version (TORAMATIC) equipped with 2 one-way valve on the connecting tube to prevent reinfusion of any aspirated fluids back into the patient.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

TORASET

SET FOR THORACENTESIS

MODEL	NEEDLE	CONNECTION	SYRINGE CAPACITY	WASTE BAG TYPOLOGY	PACKAGING
TORASET 1-55	Trio 55 mm	Stopcock 3 way	60 ml	2 lt + valve + tap	50
TORASET 1-80	Trio 80 mm	Stopcock 3 way	60 ml	2 lt + valve + tap	50
TORASET 1-VER	Veress 100 mm	Stopcock 3 way	60 ml	2 lt + valve + tap	50
TORAMATIC 1-55	Trio 55 mm	2 one-way valve	60 ml	2 lt + valve + tap	40
TORAMATIC 1-80	Trio 80 mm	2 one-way valve	60 ml	2 lt + valve + tap	40
TORAMATIC 1-VER	Veress 100 mm	2 one-way valve	60 ml	2 lt + valve + tap	40

Ordering information:

TORASET 1-55 TR 5 5 2 0 0 - 0 0TORASET 1-80 TR 8 0 2 0 0 - 0 0TORASET 1-VER TR 0 0 2 0 0 - 0 0TORAMATIC 1-55 TR 5 5 2 0 0 - 1 0TORAMATIC 1-80 TR 8 0 2 0 0 - 1 0TORAMATIC 1-VER TR 0 0 2 0 0 - 1 0

Single box: 40/50 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.