


**medax**

# VERESS NEEDLE

SPECIAL NEEDLE FOR PARACENTESIS,  
THORACENTESIS PROCEDURES  
AND PRE-LAPAROSCOPY GAS INSUFFLATION



Special needle for paracentesis, thoracentesis procedures and pre-laparoscopy gas insufflation.

- Internal retractile cannula with open or rounded tip.
- Ultra-thin needle walls to facilitate the flow of liquids.
- Penetration depth setting system.
- Penetration indicator highly visible to guarantee safe insertion.
- Available with or without stopcock and with 3-ways stopcock.
- Available with different windows type.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

# VERESS NEEDLE

SPECIAL NEEDLE FOR PARACENTESIS, THORACENTESIS PROCEDURES AND PRE-LAPAROSCOPY GAS INSUFFLATION

GAUGE	LENGTH (MM)		
14	100	120	150
15	100	120	150
16	100	120	150

Different sizes on demand

**Ordering information:**

VERESS NEEDLE WITHOUT STOPCOCK,  
OPEN TIP, INNER AND OUTER CANNULA WINDOWED **V R** GAUGE LENGTH - 0 0

VERESS NEEDLE WITH STOPCOCK,  
OPEN TIP, INNER AND OUTER CANNULA WINDOWED **V R** GAUGE LENGTH - 1 0

VERESS NEEDLE WITH 3-WAYS STOPCOCK,  
OPEN TIP, INNER AND OUTER CANNULA WINDOWED **V R** GAUGE LENGTH - 2 0

VERESS NEEDLE WITHOUT STOPCOCK,  
CLOSED TIP, 2 WINDOWS ON THE INNER CANNULA **V R** GAUGE LENGTH - 0 1

VERESS NEEDLE WITH STOPCOCK,  
CLOSED TIP, 2 WINDOWS ON THE INNER CANNULA **V R** GAUGE LENGTH - 1 1

VERESS NEEDLE WITH 3-WAYS STOPCOCK,  
CLOSED TIP, 2 WINDOWS ON THE INNER CANNULA **V R** GAUGE LENGTH - 2 1

VERESS NEEDLE WITHOUT STOPCOCK,  
CLOSED TIP, 3 WINDOWS ON THE INNER CANNULA **V R** GAUGE LENGTH - 0 2

VERESS NEEDLE WITH STOPCOCK,  
CLOSED TIP, 3 WINDOWS ON THE INNER CANNULA **V R** GAUGE LENGTH - 1 2

VERESS NEEDLE WITH 3-WAYS STOPCOCK,  
CLOSED TIP, 3 WINDOWS ON THE INNER CANNULA **V R** GAUGE LENGTH - 2 2

Single box: 20 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

