


**medax**

# VILLO-KIT

KIT FOR TRANSABDOMINAL RETRIEVAL  
OF CHORIONIC VILLI DURING PRENATAL DIAGNOSIS  
USING THE INTRODUCER CO-AXIAL NEEDLE TECHNIQUE



Kit for transabdominal retrieval of chorionic villi during prenatal diagnosis using the co-axial needle technique.

- Introducer coaxial needle for a more precise targeting during the procedure.
- Transparent Luer Lock connector.
- Sharpened cannula for a safe and painless insertion and penetration.
- Internal or external echogenic treatment for a correct positioning under ultrasound guidance.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

**medax**

# VILLO-KIT

KIT FOR TRANSABDOMINAL RETRIEVAL OF CHORIONIC VILLI DURING  
PRENATAL DIAGNOSIS USING THE INTRODUCER CO-AXIAL NEEDLE TECHNIQUE

GAUGE	COLOUR	LENGTH (MM)
20	 Yellow	20 160 20 200
21	 Light green	21 160 21 200

Ordering information:

VILLO-KIT



Single box: 20 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.